FORA G20 Test Strip

Use Only With FORA G20
Blood Glucose Monitoring System

Warnings

- For *In vitro* diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that has come into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your Blood Glucose Monitoring System Owner's Manual before you use this test strip. Use only FORA G20 Test Strips with FORA G20 Blood Glucose Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Please do not use the Glucose Monitoring System on critically ill patients. While the blood glucose result is extremely hyperglycemia (over 600 mg/dL (33.3 mmol/L)), the collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological blood glucose level might not be a true reflection. The following circumstances may apply: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemic, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA G20 test strips, when used together with FORA G20 Blood Glucose Monitoring System, allow your blood glucose levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger, and the following areas: the palm, forearm, upper arm, calf, and thigh. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

About Alternative Site Testing (AST)

IMPORTANT: There are limitations for performing AST. Please read the FORA G20 Blood Glucose Monitoring System Owner's Manual and consult your doctor before you perform AST.

Alternative site testing (AST) is when individuals check their blood glucose levels using areas of the body other than the fingertip. The FORA G20 test strips allow AST to be performed on sites other than the fingertip. We strongly recommend that you perform AST **ONLY** at the following times:

- During a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do NOT use AST if:

- You think your blood glucose is low.
- You may not notice if you are hypoglycemic.
- Your AST results are inconsistent with the way you feel.
- You are testing for hyperglycemia.
- Your routine glucose results often fluctuate.

To obtain a blood sample from the alternative sites, please **rub the puncture site for approximately 20 seconds before following** the procedures of "Testing Your Blood Glucose".

Limitations

- Hematocrit: The hematocrit level is limited to between 20% and 60%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: **This test strip must not be used for the testing of newborns.**
- Metabolites: Dopamine, L-Dopa, methyllopa, uric acid, ascorbic acid, and acetaminophen at normal blood concentration do not significantly affect blood glucose test results.
- Lipemic Effects: Blood triglycerides up to 3000 mg/dL (33.9 mmol/L) do not affect the results significantly, but may affect results at higher levels.
- Altitude Effects: Altitudes up to 10,742 feet (3,275m) do not affect test results.


Storage and Handling

IMPORTANT: Do not use the test strips if they have expired.


- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- Store the test strips in a cool, dry place between 4°C and 40°C (39.2°F and 104°F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- Store the test strips in their original vial ONLY. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Use each test strip immediately after taking it out of the vial. Close the vial immediately after taking out a strip.
- Keep the vial closed at all times.
- Do not bend, cut, or alter the test strip.

Strip Appearance

- Absorbent Hole**
Apply a drop of blood here. The blood will be automatically absorbed.
- Confirmation Window**
This is where you confirm if enough blood has been applied to the absorbent hole in the strip.
- Test Strip Handle**
Hold this part to insert the test strip into the slot.
- Contact Bars**
Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

Attention: Test results might be wrong if the contact bar is not **fully** inserted into the test slot.

Testing Your Blood Glucose

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

STEP 1
Insert the test strip fully into the slot of the meter until it will go no further. When the strip is fully inserted, the meter will do several self-checks.

STEP 2
Collect a blood sample with the test strip. A sufficient quantity of blood is required for the test to provide accurate results. Touch the blood drop with the absorbent hole of the test strip, and wait until the confirmation window is fully covered. **Do NOT** apply a smeared blood sample. The meter will start counting down.

STEP 3
After a few seconds, the meter will display your blood glucose level. The last reading will be automatically saved in the meter. Turn it off by removing the test strip and throw away the used test strip.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

Your blood glucose readings deliver plasma equivalent results and are displayed in milligrams of glucose per deciliter of blood (mg/dL).

Reference values

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meal	Less than 100 mg/dL (5.6 mmol/L)
2 hours after meals	Less than 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.

Please consult your doctor to determine a target range that works best for you.

Questionable or inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
- Make sure the confirmation window of the test strip is completely filled with blood.
 - Check the expiry date of the test strips
 - Check the performance of your meter and test strip with the control solutions.

PLEASE NOTE: Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

Questionable or inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
- Make sure the confirmation window of the test strip is completely filled with blood.
 - Check the expiry date of the test strips
 - Check the performance of your meter and test strip with the control solutions.

PLEASE NOTE: Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

Quality Control Testing

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

IMPORTANT: The reference range of the control solutions may vary with each new vial of test strips. Make sure you check the range on the label of your current vial.

Chemical Components

- > Glucose oxidase (*A. niger*) 10%
- > Electron shuttle 50%
- > Enzyme protector 8%
- > Non-reactive ingredients 32%

Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only.

Sample Size: 0.7 µL
Reaction Time: 7 seconds
System Measurement Range: 20 to 600 mg/dL
Hematocrit Range: 20% to 60%

Accuracy

The table below displays how often FORA G20 achieves this target. The chart is based on a study carried out on 118 patients to see how well FORA G20 performed compared to YSI-2300 reference method results.

For glucose results <75 mg/dL (4.2 mmol/L), the percentage (number) of FORA G20 results that were within ±15 mg/dL (0.83 mmol/L) of the reference method.	100% (19/19)
For glucose results ≥ 75 mg/dL (4.2 mmol/L), the percentage (number) of FORA G20 results were within ±20% of the reference method.	98% (97/99)

Precision

The CV (%) is less than 5% both in intermediate precision and repeatability.

Symbol Information

Symbol	Referent	Symbol	Referent
	In vitro diagnostic medical device		Batch code
	Do not reuse		Manufacturer
	Consult instructions for use		Serial number
	Temperature limitation		CE mark
	Use by		Green Dot: German dual waste disposal system
	Keep away from sunlight		Do not use if damaged
	Keep dry		Use within 6 months after opening

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